## SUITCASE® RAMP Signature Series™





## **Redesigned Top to Bottom**

Introducing the new SUITCASE® Signature Series™ ramp. While similar in appearance and function to the faithful SUITCASE Classic Series ramp, the Signature Series version features several innovative design enhancements. From the no pinch, non-protruding fulllength live hinge and durable, ergonomically designed handles to the enhanced extruded tread for superior traction, this ramp was redesigned from the ground up to be simply the best in the industry. Top Lip Extension (TLE) compatible. Available in 2' – 8' lengths with an 800-pound weight capacity.



Features a no-pinch, non-protruding live hinge for safe, silent, and maintenance free operation



Ergonomically designed flexible, nonbreakable handles offer convenient carrying and a comfortable grip



Bottom transition plates self-adjust to accommodate uneven terrain

## OPTIONAL TOP LIP EXTENSION (TLE)

The TLE will increase the reach of the SUITCASE® Signature Series™ ramp by extending the lip from the standard 3″ length to 9″. The EZ-ACCESS TLE is unique because it can be added or removed as needed, eliminating the need for a large ramp inventory. Shipped as a pair.



maintenance nee operation							
SUITCASE SS2	SUITCASE SS3	SUITCASE SS4	SUITCASE SS5	SUITCASE SS6	SUITCASE SS7	SUITCASE SS8	TLE
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
2'x 30"	3'x 30"	4' x 30"	5'x 30"	6' x 30"	7' x 30"	8' x 30"	8½" x 13 <sup>7</sup> /8" x 2 <sup>3</sup> /8"
24" x 15" x 4"	36" x 15" x 4"	48" x 15" x 4"	60" x 15" x 4"	72" x 15" x 4"	84" x 15" x 4"	96" x 15" x 4"	n/a
•	•	•	•	•	•	•	n/a
9	14	18	25	30	36	42	5 per set
800 🁛	800 🐴	800 拳	800 🀴	800 🐴	800 👛	800 🐴	**500 🐴
•	•	•	•	•	•	•	•
Curbs, 1 step	Curbs, 1 step	Curbs, 1-2 steps	Curbs, 1-2 steps, some minivans	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs	2-3 steps, some minivans & SUVs	n/a
	2'x30" 24"x15"x4" . 9 800	2'x 30" 3'x 30"  24" x 15" x 4" 36" x 15" x 4"  9 14  800 800 800					

<sup>\*</sup> Weights are approximate and may vary.

Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.





<sup>\*\*</sup> Intended for unoccupied use only.