

Be Prepared! ADA compliance required by January 2013

Confused by the new ADA regulations for fixedpool lifts? AmeriGlide makes it easier for you to comply with an attractive addition to your pool environment. AmeriGlide pool lifts are engineered for safe lifting and durability for years to come.



# AmeriGlide Pool Lift 350

# ADA COMPLIANT AMERIGLIDE POOL LIFT

## AGPOOL350 Secure Installation. Easy Removal

Concrete anchors are set 4" deep on pool deck, which easily detach for seasonal storage.

#### **Wall Mounted Battery Charger**



**Folding Foot Plate** 

#### **Durable Stainless Steel Construction**

Eliminates rust water draining into pool while also extending the life of your pool lift.

#### **Durable Finish**

Pool Lift frames are sandblasted and powder-coated to create a smooth layer of color that is then baked at temperatures up to 400° F for a long-lasting and beautiful finish. The finish will not peel, chip or fade under the sun.

### **Key Features**

- Meets and exceeds ADA specifications by lowering to 26" depth
- Fully submersible seating
- Concrete anchor mounting
- Durable stainless steel construction
- Folding arms and foot plate
- Comfortable sling seat
- Sealed hand control
- Wall-mounted charger with hanger
- Removable, easy-access battery

\*In 2010, update regulations under the Americans with Disabilities Act (ADA) adopted the 2010 Standards for Accessible Design which contain specific accessibility requirements for a number of types of recreational facilities, including swimming pools, wading pools, and spas.

In January 2012, the Department issued quidance title "ADA 2010 Revised Requirements: Accessible Pools - Accessible Means of Entry and Exit" to assist entities covered by Title III of the ADA, such as hotels and motels, health clubs, recreation centers, public country clubs, and other businesses that have swimming pools, wading pools, and spas, in understanding how the new requirements apply.

For more information on new regulations, visit www.ADA.gov .