


# Why choose the Companion?

## It's Easy to Drive

The Delta tiller and faceplate controls are intuitive and comfortable to drive. The EZ Reach tiller adjustment lever allows you to find the most comfortable tiller position.



## Co-pilot Equipped

See a flashing green light? No problem! Our convenient under-carpet troubleshooting guide will assist you through a passing storm.



## Enjoy First Class Comfort Adjustments

It doesn't get any better than this. With the plush seat, arm rests and EZ Reach tiller adjustment capabilities, you'll enjoy the most leg room and convenience as you ride or exit the scooter.



## Luxury Lighting Package

To make night time travel a breeze, we've added a new, adjustable LED headlight and ultra bright tail lights for your complete safety.



## Advanced Comfort Seating

LiquiCell® is an ultra-thin, liquid-filled pad designed to prevent pressure sores and enhance seat comfort.



GC-240

GC-340

GC-440

first class specifications

Model	GC240	GC340	GC440
Length	42.5"	45.5"	47.5"
Width	23.5"	23.5"	24"
Tires	Solid	Solid	Solid
Front Tire	9"	9"	10"
Rear Tire	9"	10"	10"
Weight Capacity	300 lbs.	350 lbs.	350 lbs.
Maximum Speed	4.5 mph	4.5 mph	4.5 mph
Ground Clearance	3.5"	3.5"	4"
Operating Range	15.5 miles*	16.5 miles*	15.0 miles*
Turning Radius	39"	40.5"	53.5"
Batteries (2)	Type U-1	Type U-1	Type U-1
Total Weight w/ Batteries	166 lbs.	172 lbs.	188 lbs.
Standard Seat Stadium Style	18" x 16"	18" x 16"	18" x 16"
Freewheel Mode	Yes	Yes	Yes
Electronic Speed Control	Yes	Yes	Yes
Brakes	Electromechanical	Electromechanical	Electromechanical
Charger	Off-board	Off-board	Off-board
Rear-wheel Drive	Sealed transaxle motor, 24 VDC	Sealed transaxle motor, 24VDC	Sealed transaxle motor 24VDC
Controller	90 amp	90 amp	90 amp
Warranty- Frame	Lifetime Limited	Lifetime Limited	Lifetime Limited
Warranty - Motor/Electronics**	2 Years	2 years	2 years
Warranty - Drive Train	2 Years	2 years	2 years

## optional accessory upgrades



Comfort Ride



Cup Holder



Safety Flag



Rear Basket



Crutch Holder



Power Elevating Seat



Oxygen Holder



Walker Holder



Pack'n'Go



Quad Cane Holder

# TRAVEL FIRST CLASS



under-carpet troubleshooting guide

LiquiCell® Comfort Equipped seat

auto lock tiller pin

seat slide adjustment lever

EZ Reach tiller adjustment lever

Delta tiller

easy to read control panel

enhanced tiller design

foldable stadium seat

black, non-marking tires

removable shopping basket with handle

LOTS of leg room

4 mobility transportation anchors

low profile tires on high-strength metal rims

charging port

auto lock tiller pin

optional auto ride comfort seat

width & height adjustable armrests

REALLY COOL WHEELS!

comfortable padded armrests

excellent turning ability

seat rotation lever

adjustable LED headlight

washable automobile-style carpeting

comfortable padded armrests

THE GOLDEN COMPANION IS YOUR TICKET TO RIDE FIRST CLASS EVERY DAY!

www.goldentech.com

**GOLDEN**  
technologies  
Your Life in Motion

800.624.6374

INTRODUCING THE MOST ADVANCED SCOOTER IN ITS CLASS

**GOLDEN**  
**COMPANION**



## KEEP CONNECTIONS WITH THOSE YOU LOVE

Allow the GC-240 to enhance your life and take you the places you've always had trouble getting to. This midsize scooter is ideal for trips to the mall, dining out or an afternoon at the park with family.

### STANDARD FEATURES ON ALL COMPANION MODELS INCLUDE:

- Delta tiller for easy, comfortable driving
- LiquiCell® comfort-equipped seat
- Stylish, low profile tires on strong metal rims
- Width, height, and angle adjustable armrests



**GC-240**



## REACH HIGHER ALTITUDES!

The GC-340, shown with the optional power elevating seat, gives you more legroom and plenty of style to be on the go indoors and out. Its impressive leg room allows you to stretch your legs and ride in comfort even on the longest tours of your favorite destinations.

### STANDARD FEATURES ON ALL COMPANION MODELS INCLUDE:

- Removable shopping basket with handle
- Side view mirrors for safe navigation
- Easy to reach charger port
- Sliding seat design for 4" of adjustment



**GC-340**



## MAINTAIN YOUR CRUISING SPEED AS YOU NAVIGATE THE GREAT OUTDOORS

Experience the great outdoors like you've always wanted to with the four-wheel GC-440. With 10" tires on the front and rear, the GC-440 has an impressive 4" ground clearance and the durability to traverse over many obstacles.

### STANDARD FEATURES ON ALL COMPANION MODELS INCLUDE:

- Angle-adjustable LED headlight for great visibility
- Washable, automobile-style carpeting
- Convenient under-carpet troubleshooting guide
- The new "EZ Reach" tiller adjustment lever



**GC-440**

