

F U L L S I Z E S C O O T E R



Seat Color
Black Vinyl

Body Colors
Red

Blue

GOLDEN COMPANION

3 wheel
GC340C

The Companion three-wheel scooter offers more leg room and foot room with a wider front deck and a stylish new design. Features include easy disassembly; an adjustable LED headlight; ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion is the ultimate value in a full-time portable scooter.



GC-340C

Companion 3 wheel Specifications

Model	GC340C
Medicare Code	Pending
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5.3 mph
Operating Range ¹	20.9 miles (U1) 25.4 miles (40AH)
Ground Clearance Under Deck	4.5"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	41"
Type Batteries	2-U1, (2-40AH optional)
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
Dimensions:	
Length	45.5"
Width	23.5"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	24" - 26"
Ground to Top of Deck	6.25"
Deck to Top of Seat	18" - 20"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.5"
Weight of Unit: (assembled)	176 lbs (U1) 192 lbs (40AH)
Front Half	56 lbs
Rear Half	48 lbs
Seat w/Arms	26 lbs
Batteries	21 lbs ea. (U1), 29 lbs ea. (40AH)
Tires:	Solid/Low Profile
Front	9"
Rear	10.75"
Rear Anti-Tippers	3"
Standard Seat:	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" X 16"
Color	Black Vinyl
Optional Seats:	Stadium Style
Back Height - no headrest	16"
Width x Depth	20" x 18"
Color	Black Vinyl
Optional Seats 2:	Captain's Seat
Back Height - no headrest	19"
Back Height with headrest	29"
Width x Depth	18" x 18"
Color	Black Vinyl
Seat Options:	
Pneumatic Rise Seat	Yes
Power Elevating Seat	Yes
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	1 Year



Angle adjustable LED headlight creates a wide path of bright light



Quiet & maintenance-free motor




EZ Reach tiller adjustment



Patented infinite adjustable tiller



Disassembles easily!



Standard Features on Every Model:

- Under mat troubleshooting guide
- Adjustable sliding seat
- Delta tiller with auto lock mechanism
- Adjustable LED headlight
- Easy reach tiller adjustment
- U-1 batteries standard

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.
² Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

