



PR458 shown in
Buckskin



Admiral



Palomino



Buckskin



Coffee Bean

PR458

If you are looking for a simple power lift recliner with a modern twist, then look no further. This three position recliner features a stylishly sleek arm, seam back, and full chaise seat for ultimate comfort and support. In addition to its fashionable and relaxing design, this chair will lift you from a seated to a standing position with the touch of a button.

PR458

Model Number	PR458-SMA	PR458-MED	PR458-LAR	PR458-S23	PR458-M23
Size	Small	Medium	Large	Small 23" wide	Medium 23" wide
Weight Capacity	300 LBS	375 LBS	375 LBS	375 LBS	375 LBS
Recline Positions	3	3	3	3	3
Standard Fabrics	4	4	4	4	4
Fabric Options	None	None	None	None	None
Chaise Seat	YES	YES	YES	YES	YES
Coil Springs	NO	NO	NO	NO	NO
Back Type	Button	Button	Button	Button	Button
Head Rest Cover	YES	YES	YES	YES	YES
Arm Covers	YES	YES	YES	YES	YES
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	30"	32"	34"	34"	34"
Overall Height	40"	43"	44"	40"	43"
Width Between Arms	19"	20.5"	23.5"	23"	23.5"
Floor to Top of Seat	18"	20"	21"	18"	20"
Seat Depth	19"	20"	22"	19"	20"
Seat to Top of Back	27"	28"	29.5"	27"	28"
Seat to Top of Arm	7.5"	6.5"	7.5"	7.5"	6.5"
Seat Angle Raised	15°	15°	15°	15°	15°
Floor to Top of Foot Rest Reclined	21.5"	23"	23"	21.5"	23"
Floor to Top of Back Reclined	32.5"	33.5"	34"	32.5"	33.5"
Back Angle Closed	107°	107°	107°	107°	107°
Back Angle Reclined	145°	145°	145°	145°	145°
Foot Rest Angle Closed	90°	90°	90°	90°	90°
Foot Rest Angle Opened	175°	175°	175°	175°	175°
Distance Required from Wall Reclined	17"	17"	18"	17"	17"
Extension reclined from Back of Seat to end of Foot Rest	40.5"	42.5"	43.5"	40.5"	42.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	NO	NO	NO	NO	NO
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.	1 Yr.	1 Yr.
Removable Back Option Available	NO	NO	NO	NO	NO
Smart Tek System Standard	YES	YES	YES	YES	YES
Customizable	NO	NO	NO	NO	NO
Left Hand Control Option Available	NO	NO	NO	NO	NO
Extra Pocket Option Available	NO	NO	NO	NO	NO
Heat & Massage Option Available	NO	NO	NO	NO	NO
Gel Memory Foam Option Available	NO	NO	NO	NO	NO
Power Pillow Option Available	NO	NO	NO	NO	NO
AutoFlex Option Available	NO	NO	NO	NO	NO
Foot Rest Extension Option Available	NO	NO	NO	NO	NO
Straight Lift Option Available	NO	NO	NO	NO	NO
User Height (suggested)	5'4" - 5'10"	5'4" - 5'10"	5'11" - 6'2"	5'4" - 5'10"	5'4" - 5'10"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Literature is current at time of printing. AmeriGlide reserves the right to make changes to the product or literature at any time.