AmeriGlide 325M

3 POSITION LIFT CHAIR

- 325 lb. Weight Capacity
- Full Chaise Pad
- Battery Backup
- Button Back
- Convenient Remote Control
- Made in America

There is no need to settle for limited mobility. The AmeriGlide 325M provides a stylish chaise seat with a classic design. Choose the color option that best suits your home and you can relax while also improving your well-being.

AmeriGlide Accessibility Solutions



MADE IN THE USA

Literature is current at time of this version. AmeriGlide reserves the right to make changes to the product or literature at any time.

325M LIFT CHAIR

The AmeriGlide 325M 3 Position lift chair offers a weight capacity of 325 pounds and a three position recline capability. The 325M 3 Position features a full chaise pad, is backed by a one-year in-home service warranty, and is available in an array of colors to match your home decor. The 325M is built with an integrated battery backup to keep you mobile even if your home loses power, and its hardwood frame ensures that your chair will remain sturdy and hold up even after repeated use.



CAMEL

FEATURES

- · Fabrics Available: Camel, Chocolate, Navy & Charcoal
- 1 Year In-Home Parts and Labor
- Battery Backup

SPECIFICATIONS

Recline Positions	3
Weight Capacity	325 lb
Min. Recommended Height	5' 4"
Max. Recommended Height	5' 8"
Overall Height	42.5"
Overall Width	33"
Seat Width	19"
Distance Required from Wall Reclined	18"
Seat to Top of Back	25"
Seat Depth	19"
Floor to Top of Seat	18"

NOTE: Measurements may vary +/- 1". Specifications subject to change without prior notice.

Literature is current at time of this version. AmeriGlide reserves the right to make changes to the product or literature at any time.





CHOCOLATE



NAVY