

Trouble lifting and transporting a heavy power chair? Now you can!



High Capacity 2-Axis Hoist Lift

Fully-powered lifting and rotation for effortless loading of a 400 lb+ power chair.

Features & Benefits

- Ideal for pick up trucks and full size vans
- Heavy duty construction and stronger components to lift fully-equipped power wheelchairs
- Weather-resistant hand control included; optional wireless remote control available
- Load/unload mobility device conveniently from side of vehicle
- Passenger or driver-side install (DS02 option)
- Adjustable base adapters for faster installation

Specifications	AL425HD Axis II Lift
Accommodates	Heavy power chairs
Lifting Capacity	450 lb (204 kg)
Installed Weight	70 lb (32 kg)
Adjustable Boom Height	31.5" – 46.5"
Adjustable Boom Length	26.25" – 33.5"
Adjustable Boom Pitch	0° - 15°
Securement	Docking device included
Options	Battery pack, wireless remote, folding arm
Warranty	3-year transferable