

MODEL

Jazzy Select 6-SS ¹	Group 2: Standard (K0822)
Jazzy Select 6-SS ¹	Group 2: Standard (K0823)
Weight Capacity	300 lbs.
Turning Radius	21.875"
Overall Width ³	23.5"
Overall Length ³	35.5"
Maximum Speed ²	Up to 4 mph
Ground Clearance ³	1.75"
Front Wheels	6" solid casters
Drive Wheels	10" solid
Rear Wheels	6" solid casters
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Active-Trac® (ATX) with Mid-wheel 6® Technology
High-back Seat Sizes	16"x16"-18", 18"x18"-20", 20"x18"-20"
Seat-to-Floor Range	21.25" - 23.25"
Per Charge Range ^{3,5}	Up to 15 miles
Standard Electronics	40A, PG GC3*
Component Weights	Base 97.4 lbs. (<i>heaviest piece</i>) Seat 50.2 lbs. Batteries 25 lbs. each
Battery Requirements ^{4,6}	(2) 12 volt, U-1 batteries
Battery Charger	3.5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

¹ Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

⁴ AGM or gel cell type required.

⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁶ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

*Required to upgrade the controller to VR2 6-key controller.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

0002615-021521

Jazzy^{® 1}
SELECT 6

**SPORTY DESIGN
AND SUPERIOR
PERFORMANCE**



FDA Class II
Medical Device¹

AMERICA'S FAVORITE POWER CHAIR

gojazzy.com |    

A UNIQUE RIDE EXPERIENCE



Depth-adjustable, limited-recline, comfort high-back seat

Height, width and angle adjustable armrests

Intuitive controller is easy to operate

Flip-up, height-adjustable foot platform offering two positions of height and one position of depth

Color-through shroud



Active-Trac® Suspension with Mid-wheel 6® Technology




PG GC3 controller

The **Jazzy Select® 6¹** incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac® Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions.

- Mid-wheel drive technology and a compact size is perfect for small spaces by offering a tight turning radius
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Easy-access free-wheel lever
- Color-through shroud keeps the shroud looking new by reducing the visible scrapes and chips

Available Colors



Jazzy Red Jazzy Blue

Seating

Black, vinyl, depth-adjustable and limited recline standard comfort high-back with 60" lap belt available sizes:

- 16"W x 16" - 18"D
- 18"W x 18" - 20"D
- 20"W x 18" - 20"D

Optional Seating

- Black vinyl, high-back seat up to 20" x 20"

Options

- Bell-shaped foot platform
- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away footrests
- Swing-away joystick
- VR2 6-key (must and only upgrade with PES option)

Accessories

- Cane/crutch holder
- Cup holder
- Flag
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover—medium
- XLR USB