SUMMIT SL3500D

OUTDOOR CONVENIENCE







Large weatherproof cover included

Folds to 14" from wall

The SL350OD outdoor stair lift is engineered with a durable extruded aluminum track and marine-grade vinyl seating to handle rain, wind, snow, and temperatures from 0°F to 125°F. Includes header, internal chassis, and a large weatherproof cover for protection from outside elements.

FEATURES AND BENEFITS

- Designed for long-term outdoor use
- Header cover, internal chassis cover and a large weatherproof cover included
- Rocker controller in armrest operates with constant-pressure functionality to stop the lift immediately if the control is released
- Extruded aluminum rail mounts to the stairs (not the wall) on either side of staircase
- Operates on standard 115V AC household current; no batteries required
- Can be installed at 350 lb capacity up to 15'9" and 250 lb up to 25'
- Call/send controls at top and bottom of stairs





My wife had trouble walking up to our front door from the street. Having the SL350OD allows us to enjoy our forever home for years to come by giving her access to the home year around.

SL350OD PRODUCT SPECIFICATIONS

WEIGHT CAPACITY 350 lb up to 15' 9" track 250 lb for 15'10" - 25' track

TRACK LENGTH 16' to 25'

SEAT SIZE 17.5" x 14"

FOLDED WIDTH 14" from wall

- Footrest obstruction switch

- Swivel seat cut-off switch

SAFETY - Final limits

FEATURES - Slack cable brake with switch

- Constant pressure controls

- Seat belt

INCLINE LIMITS 30° to 45°

SPEED 20 FPM

FOOTREST SIZE 9.5" x 10"

DRIVE Cable

WARRANTY 3-year parts

OPTIONS

Key lock

Extra rail



