

GP302 GOLDEN CRICKET

The all-new Cricket Power Wheelchair is stylish, agile, and extremely portable. Crafted with an incredibly strong yet lightweight 100% carbon fiber frame, the Cricket weighs as light as 33 lbs.*, making travel a breeze! The slim sporty design is perfect for navigating indoors and out while sightseeing, shopping, or dining. The roomy padded seat cushion and durable mesh back align with the natural curves of your body providing long-lasting comfort while on the go. The Cricket features a well-styled joystick designed to make everyday tasks easier such as driving up to a table or maneuvering through a doorway. It has an impressive 22" turning radius and offers a 300 lb. weight capacity, maximum speed of 3.7 mph, and an operating range of up to 9.3 miles on airline-approved 12 AH 24V lithium-ion batteries. It also offers convenient offboard charging and folds effortlessly for transport. The Cricket offers excellent performance and portability for your active lifestyle!



Standard Color



Carbon Fiber



Standard Features:

- 100% Carbon Fiber Frame
- 300 lb. Weight Capacity
- Removable Compact
 Joystick with USB Charger
- 35 Amp Controller
- Convenient Under Seat Storage Bag
- Removable Padded Seat
- Offboard Charging Capability
- Easy One-Step Folding
- Stands on End for Space-saving Storage
- Airline-approved Lithium-ion Battery



Model Number	GP302
Medicare Code	K0899 Not Coded
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	3.7 mph
Operating Range ¹	9.3 miles
Ground Clearance	2.6"
Turning Radius	22"
Type Batteries	12AH 24V Lithium-Ion
Color	Black Carbon Fiber
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	35 Amp
Dimensions:	
Length	36.6"
Width	22.25″
Height (ground to top of back of seat)	36.5″
Ground to Top of Seat	19.3″
Number of Seat Height Positions	N/A
Size of increments	N/A
Weight of Unit: (assembled)	39.3 lbs.
Weight of Unit: (minus parts)*	33.2 lbs.
Front Half	N/A
Rear Half	N/A
Seat w/Arms	N/A
Batteries	3.8 lbs.
Tires:	Flat Free
Drive Wheels (Rear)	8.5″
Casters (Front)	7"
Rear Anti-Tip Wheels	1.6″
Standard Seat:	
Back Height - no headrest	19.3″
Width x Depth	17.3" x 16.5"
Color	Black
Footrest:	Yes
Height adjustable	N/A
Angle adjustable	N/A
Warranty:	
Frame	5 Years
Drive Train	13 Months
Electronics ²	13 Months
Batteries	12 Months



Sporty, strong and high-performing with a 100% carbon fiber frame!

FOLD & STOW! Compact Dimensions 30.2" × 22.25" × 14"





A T

Convenient USB Flip Up

Charging Port









Removable 1.25" Thick

Padded Seat

Stands on End for Storage



Phone: 800-624-6374 Fax: 800-628-5165 www.goldentech.com

Find us on Facebook www.facebook.com/goldentechnologies

 Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.
 ² Electronics warranty excludes batteries.

* Excludes battery, joystick, cushion and bag.

FDA Class II Medical Device

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

