

VICTORY[®]

PLATINUM

PROUDLY MADE IN
AMERICA
OF FOREIGN AND DOMESTIC PARTS

EXPERIENCE THE PLATINUM DIFFERENCE



FDA Class II
Medical Device¹



FDA Class II
Medical Device¹



PRIDE
MOBILITY[®] pridemobility.com

1) The Pride[®] Victory Platinum is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

VICTORY® PLATINUM


SC710PLA



FDA Class II Medical Device¹

FEATURES

- **MADE IN AMERICA** from foreign and domestic parts
- Patented **EZ TURN** technology provides an impressive turning radius of only 45.5"
- **EXCLUSIVE CHARGE360** magnetic charger
- Easy one-handed tiller adjustment
- **FULL FRONT AND REAR SUSPENSION** with optional suspension seat post
- Feather-touch disassembly for convenient transport or storage
- 18x17" pinchless hinge, black vinyl seat with 8-point swivel & sliders
- Maximum intensity LED headlight
- Front basket and pair of rear view mirrors (standard)
- Water-resistant console



Exclusive **Charge360** magnetic charger



Easy one-handed tiller adjustment



Patented **EZ Turn** technology provides unrivaled tight turning.

ACCESSORIES

(Only one rear accessory can be used at a time)

- Cane Clip Attachment
- Cane Crutch Holder
- Cell Phone Holder
- Cup Holder
- Forearm Crutch Holder
- Lap Belt
- Military Patches
- Oxygen Tank Holder
- Rear Basket
- Saddlebag
- Safety Flag
- Single/Double Cane Holder
- Suspension seat post
- Walker Holder
- Weather Cover
- Wishbone Crutch Holder
- XLR USB Charger

- *Cup holder, rear basket, and saddlebag are available as a bundle.*
- *High back seats with limited recline are available as upgrade options.*

SPECIFICATIONS

WEIGHT CAPACITY	400 lbs.
MAXIMUM SPEED ²	Up to 6 mph
DIMENSIONS ³	Length: 47" Width: 24"
TURNING RADIUS	45.5"
RANGE ^{3,4}	Up to 15.5 miles

WARRANTY

FRAME:	Lifetime limited
ELECTRONICS:	2-year limited
DRIVE MOTORS:	2-year limited

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments. K0899 Power Mobility Device not coded by DME PDAC or does not meet criteria. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

MATTE

Iceberg Blue

MATTE

Cardinal Red

AUTHORIZED PROVIDER



pridemobility.com

PRIDE
MOBILITY®