Live fully independent with the **Infinity Curved Stair Lift**

Stay living independently in your home.

Your stairs should never be your biggest obstacle while aging in place. Negotiating the stairs shouldn't stop you from doing the things you want to do.

By choosing the Infinity curved stair lift, you can continue enjoying your entire home without worry.

Give peace of mind to yourself and others by installing an Infinity in your home.

With minimum assistance from others and its made-to-measure rail, the Infinity stair lift will keep you in good hands.

A Rail Designed Around Your Home

The Infinity rail is specifically designed to fit your staircase, no matter the shape.

The durable powder coated rail gives stability and ensures a smooth ride around each corner giving the best ride possible.



Color Options

The Infinity has a variety of vinyl seat covers to choose from to match your style and home.



087_InfinityCurvedSL_OneSheet_rev122723

Technical Specs

	mm	inch	DF
A Overall Height (from floor)	1058	42	
B Height of arms (from floor)	837	33	
C Height of seat (from floor)	594	23	
D Overall width	607	24	
E Width of footrest	300	12	
F Overall depth (from wall)	620	25	
G Depth of seat (to backrest)	385	15	
H Depth of footrest	270	11	
Overall Height (when folded)	1078	43	
J Overall depth (when folded)	407	16	

Standard Features

The Infinity stair lift include the following as standard.



Retractable Seatbelt

Just like a car seatbelt, easy to adjust and use, giving the user the best comfort and safety possible whilst operating the stair lift.



Locking Keys

Two keys are included for locking the stairlift in an off position for added safety, preventing unwanted operation.



Remote Controlled

Two remote controls are supplied giving convenient operation. The stair lift can be called up the staircase when it's not in used.



Comfort and Control

The soft start/stop function is enhanced by the complete comfort of the padded seat, backrest and ergonomic controls.

Why Choose the Infinity curved?

- One of the slimmest rails
- Smooth Start / Stop Function
- Fold Away Design
- Power Swivel

- Removable Seat Covers for Easy Cleaning
- Stop Motion Detector
- Easy Swivel Lever Controls
- Ergonomic Controls

